

## Bone Graft Substitutes

# MIS System



### Control, Precision, Reliability.

The MIS System represents an industry-first cartridge delivery system, catering to the minimally invasive surgery field. This innovative system allows for delivering pre-loaded cartridges of NovaBone Putty, the next-generation synthetic bone graft.

### Properties:

- MIS System ensures controlled, precise, and reliable graft delivery
- Preloaded NovaBone putty cartridges with volume markings
- Requires no preparation and can be directly applied to the intended graft site
- Flexible tip allows for graft placement in difficult-to-reach areas

### Indications:

- Indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure
- Indicated to be gently packed into bony voids or gaps of the skeletal system. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone

### MIS System

	USA	International
• MIS Putty System, 10cc .....	NB6610 .....	772610
• MIS Putty Refill Cartridge, 5cc .....	NB6650 .....	772605
• MIS Dispenser .....	NB6600 .....	772600
• MIS Putty System, 15cc .....	N/A .....	772615

