

Bone Graft Substitutes

NovaBone Putty®



Saving Time. Enhancing Precision.

NovaBone Putty® is a versatile bone graft substitute renowned for its exceptional handling characteristics, saving valuable time and enhancing placement precision. Its specially formulated consistency allows for easy packing into voids of any size or shape without sticking to surgical gloves.

Properties:

- Ready to use, requires no preparation, and can be directly applied to the intended graft site
- Non-hardening putty with exceptional handling characteristics
- Resists migration during irrigation
- Does not adhere to surgical gloves
- Highly malleable, easy to mold, and pack or graft directly into any defect
- Provides a bone void filler that resorbs and is replaced with bone during the healing process

Indications:

- Indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure
- Indicated to be gently packed into bony voids or gaps in the skeletal system (i.e. the extremities, posterolateral spine, pelvis, and intervertebral disc space). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone.

NovaBone Putty®

	USA	International
• 0.5cc Novabone Putty CMF	N/A	770600*
• 1.0cc Novabone Putty CM.....	N/A	770601*
• Putty Clamshell 2.5cc	NB0602	770602*
• Putty Clamshell 5.0cc	NB0605	770605*
• Putty Clamshell 10cc	NB0610	770610*
• Putty Syringe 0.5cc	N/A	771600*
• Putty Syringe 1.0cc	NB2601	771601*
• Putty Syringe 2.5cc	NB2602	771602*
• Putty Syringe 3.0cc	N/A	771603*
• Putty Syringe 5.0cc	NB2605	771605*
• Putty Syringe 10cc	NB2610	771610*
• Putty Tapered Syringe 2.5c.....	NB2702	772702
• Putty Tapered Syringe 5.0cc	NB2705	772705



*This product is available in the EU