

NOVABONE[®]

BIOACTIVE SYNTHETIC GRAFT

MIS Cartridge Delivery System



Controlled, Exact, Reliable Delivery

NovaBone MIS delivers the only synthetic graft that signals and recruits osteoprogenitor cells while controlling the cell cycle to favor proliferation and differentiation of cells that generate new bone. NovaBone demonstrated equivalent rates of bone growth when compared to autograft.



NOVABONE[®]

BIOACTIVE SYNTHETIC GRAFT

MIS Cartridge Delivery System

- The NovaBone Putty MIS Cartridge Delivery System comes preloaded, and is available with multiple cartridges for ease of use.
- Flexible tip allows for graft placement in difficult to reach areas.
- Discrete volume dispensing, 1/4cc per handle cycle.
- Plunger marking for ease of monitoring dispensed volume.
- Approved for use as stand-alone device in posterolateral spine.



MIS Cartridge Delivery System

Reference No.	Description
NB6610	MIS Cartridge System, 10cc
NB6600	MIS Cartridge Handle
NB6650	MIS Cartridge Refill, 5cc

Customer Service 855-265-8013

13631 Progress Blvd., #600
Alachua, FL 32615
386-462-7660
Facsimile 386-418-1636